

Practice Problem #15  
PHY 540 - Spring 2018

Use the following procedure to construct a VI that uses the ADR112 serial board to create a voltmeter to read single-ended 12-bit voltages in the range from 0 to 5V:

1. Use the FP and BD given here as a guide in constructing your VI.
2. Connect the ADR112 serial board to the computer with the USB-to-Serial converter cable.
3. Use one frame of a Flat Sequence Structure on the BD to enclose your voltmeter VI.
4. Construct your VI to take voltage readings every 250 ms with a While Loop on the BD.
5. Display the voltages on the FP with a Numeric Indicator and a Meter.
6. Change the scale on the Meter so it represents voltage readings from 0 to 5 V.
7. Set the Numeric Indicator on the FP to read to 3 decimal places.
8. Have an LED on the FP that will light as a warning if the voltage gets greater than 4.5 V.
9. Have a control on the FP to select the computer COM Port.
10. Have a control on the FP to select analog input Channel 0 or 1 on the ADR112 board.
11. Have a stop switch on the FP.
12. Document and create an icon for your VI.
13. Print a copy and save an electronic copy of your VI for your notebook.

Some components you need for this VI are at the following locations:

Instrument I/O >> Serial >> VISA Configure Serial Port  
Instrument I/O >> Serial >> VISA Write  
Instrument I/O >> Serial >> VISA Read  
Instrument I/O >> Serial >> VISA Close  
Programming >> Array >> Index Array  
Programming >> Strings >> Carriage Return Constant  
Programming >> String >> Concatenate Strings  
Programming >> Dialog & User Interface >> Simple Error Handler.vi  
Programming >> String >> String/Number Conversion > Decimal String to Number  
Programming >> Cluster, Class & Variant >> Unbundle By Name  
Programming >> Numeric >> Expression Node  
Programming >> Numeric >> Conversion >> To Double Precision Float

